

VIRYA BANDAGE SUPPLIES



Supplier of Surgical Gauze & Bandage

Customer Care- India: +91 85250 70186
Email: info@viryabandagesupplies.com

Australia: +61 488 499 556
www.viryabandagesupplies.net



Products We Offer

Absorbent Gauze & Bandage Cloths

Cotton Roller Bandage

Dressing Pads (Sterile & Non Sterile)

Gauze Swabs (Sterile & Non Sterile)

Abdominal Sponges (Sterile & Non Sterile)

Gamjee Rolls

Crepe Bandages

POP Bandages

Soft Rolls (Orthopedic Pad Casting)

Surgical Cotton Rolls

Surgical Cotton Balls

Disposable Bed Sheets with Pillow Covers

Under pads

Othopedic products

ABSORBENT GAUZE

DESCRIPTIONS:

Absorbent gauze has a loose, open weave, which allows fluids from the wound to be absorbed into the fibers, wicked away, or passed through into other absorbent materials in the wound's dressing.



Usage Applications	Hospital & Clinical
Color	White
Feature	Wound Dressings
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Quality	BP Typepe 17 & 13, Super Fine (42 TPI,I, 24 x 18 19 x 15 Mesh, 40s x 40s) SCH F(II)

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	90 cm x 10 mtr	100 PKTS
2	100 cm x 10 mtr	100 PKTS
3	90 cm x 16 mtr	100 PKTS
4	100 cm x 16 mtr	100 PKTS
5	90 cm x 18 mtr	100 PKTS
6	100 cm x 18 mtr	100 PKTS
7	90 cm x 20 mtr	100 PKTS
8	100 cm x 20 mtr	100 PKTS
9	120 cm x 10 mtr	100 PKTS
10	120 cm x 15 mtr	100 PKTS

ABSORBENT BANDAGE CLOTH



DESCRIPTIONS:

A bandage is a piece of fabric or other material that is used to cover, protect, or support a wound or injured body part. Bandages are typically made of materials such as cotton, gauze, or elastic, and may be adhesive or non-adhesive. They are often used in conjunction with other medical treatments, such as wound cleaning or the application of medication.

Usage Applications	Hospital & Clinical
Color	White
Feature	Wound Dressings
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Quality	Premium (60 TPI 38 x 22 Mesh, 20s x 30s, 40s x 40s) SCH F (II)

AVAILABLE SIZES :


(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	100 cm x 10 mtr	100 PKTS
2	100 cm x 20 mtr	100 PKTS

COTTON BANDAGE (ROLLER BANDAGE)

DESCRIPTIONS:

A roller bandage is a type of bandage that comes in the form of a long, continuous strip of material wound onto a roll or spool. Roller bandages are commonly made of materials such as cotton or elastic, and may be adhesive or non-adhesive. They are often used to wrap around a body part or injury site, such as an arm or leg, to provide support, compression, or protection.




Usage Applications	Hospital & Clinical
Color	White
Feature	Wound Dressings
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Quality	BP Type 13 & 17, Fine (34 TPI 19 x 15 & 24 x 18 Mesh, 40s x 40s) SCH F(II)

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	5 cm x 3 mtr	300 PKTS
2	7.5 cm x 3 mtr	250 PKTS
3	10 cm x 3 mtr	200 PKTS
4	15 cm x 3 mtr	100 PKTS

5	5 cm x 4 mtr	300 PKTS
6	7.5 cm x 4 mtr	250 PKTS
7	10 cm x 4 mtr	200 PKTS
8	15 cm x 4 mtr	100 PKTS



<i>S.No</i>	<i>SIZE</i>	<i>MOQ</i>
9	5 cm x 8 mtr	100 PKTS
10	7.5 cm x 8 mtr	75 PKTS
11	10 cm x 8 mtr	50 PKTS
12	15 cm x 8 mtr	50 PKTS

13	5 cm x 10 mtr	100 PKTS
14	7.5 cm x 10 mtr	75 PKTS
15	10 cm x 10 mtr	50 PKTS
16	15 cm x 10 mtr	50 PKTS

17	5 cm x 12 mtr	50 PKTS
18	7.5 cm x 12 mtr	50 PKTS
19	10 cm x 12 mtr	50 PKTS
20	15 cm x 12 mtr	50 PKTS

STERILE GAUZE SWABS



DESCRIPTIONS:

sterile gauze swabs are widely used to clean and cover minor wounds, to absorb secretion and to treat secondary healing wounds. Can be used to clean or cover minor wounds, to absorb light exudates.

Usage Applications	Hospital & Clinical
Color	White
Feature	Dressing
Material	100% Cotton, Soft & Hygienic
Sterility	Sterile ETO
Quality	BP Type 17, Super Fine (42 TPI 24 x 18 Mesh, 40s x 40s Count)

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	5 cm x 5cm x 8 ply	450 POUCH
2	7.5 cm x 7.5 cm x 8 ply	450 POUCH
3	10 cm x 10 cm x 8 ply	450 POUCH
4	5 cm x 5cm x 12 ply	500 POUCH
5	7.5 cm x 7.5 cm x 12 ply	500 POUCH
6	10 cm x 10 cm x 12 ply	500 POUCH

NON STERILE GAUZE SWABS



DESCRIPTIONS:

sterile gauze swabs are widely used to clean and cover minor wounds, to absorb secretion and to treat secondary healing wounds. Can be used to clean or cover minor wounds, to absorb light exudates.

Usage Applications	Hospital & Clinical
Color	White
Feature	Dressing
Material	100% Cotton, Soft & Hygienic
Sterility	Sterile ETO
Quality	BP Type 17, Super Fine (42 TPI 24 x 18 Mesh, 40s x 40s Count)

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	5 cm x 5cm x 8 ply	300 POUCH
2	7.5 cm x 7.5 cm x 8 ply	200 POUCH
3	10 cm x 10 cm x 8 ply	100 POUCH
4	5 cm x 5cm x 12 ply	200 POUCH
5	7.5 cm x 7.5 cm x 12 ply	150 POUCH
6	10 cm x 10 cm x 12 ply	100 POUCH

STERILE ABDOMINAL GAUZE SPONGES



DESCRIPTIONS:

They are used mostly during abdominal and thoracic surgery or in deep wounds to absorb fluids; however, they can also be used to “wall off” the surgical site. The friction of their coarse fibers also allows laparotomy sponges to be used to retract slippery visceral organs

Usage Applications	Hospital & Clinical
Color	White
Feature	Surgical Abdominal Dressing
Material	100% Cotton, Soft & Hygienic
Sterility	Sterile ETO
Quality	BP Type 17, Super Fine (42 TPI 24 x 18 Mesh, 40s x 40s Count)

AVAILABLE SIZES :


(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	25 cm x 25 cm x 8 ply	300 POUCH
2	25 cm x 40 cm x 8 ply	300 POUCH
3	30 cm x 30 cm x 8 ply	300 POUCH

NON STERILE ABDOMINAL GAUZE SPONGES

DESCRIPTIONS:

They are used mostly during abdominal and thoracic surgery or in deep wounds to absorb fluids; however, they can also be used to “wall off” the surgical site. The friction of their coarse fibers also allows laparotomy sponges to be used to retract slippery visceral or- gans



Usage Applications	Hospital & Clinical
Color	White
Feature	Surgical Abdominal Dressing
Material	100% Cotton, Soft & Hygienic
X-Ray	Yes
Tail	Yes
Sterility	Non - Sterile
Quality	BP Type 17, Super Fine (42 TPI 24 x 18 Mesh, 40s x 40s Count)

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	25 cm x 25 cm x 8 ply	100 PKT
2	25 cm x 40 cm x 8 ply	100 PKT
3	30 cm x 30 cm x 8 ply	100 PKT
4	25 cm x 25 cm x 12 ply	100 PKT
5	25 cm x 40 cm x 12 ply	100 PKT
6	30 cm x 30 cm x 12 ply	100 PKT

GAMJEE ROLL



DESCRIPTIONS:

Gamjee roll are highly absorbent cotton are available in low adherent gauze cover and is applied over dressing to provide low adherent wound interface. Generally used for primary wound dressings or is applied over a dressing, these bandage roll allow low adherent wound interface.

Usage Applications	Hospital & Clinical
Color	White
Feature	Dressing
Material	100% Cotton, Soft & Hygienic
X-Ray	Non Sterile
Quality	BP Type 17, Super Fine (42 TPI 24 x 18 Mesh, 40s x 40s Count)

AVAILABLE SIZES:



(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	10 cm x 3 mtr	100 ROLLS
2	15 cm x 3 mtr	100 ROLLS

GREPE BANDAGE ROLLS

DESCRIPTIONS:

Crepe bandage is suitable for use as support for Sprains and strains in joints and muscles. It also Works as a moderate compression bandage for Knee swelling, ankle swelling, and other relevant Injuries



Usage Applications	Hospital & Clinical
Color	Skin
Feature	Ortho Care
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Brand	Osseo Care

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	6 cm x 4 mtr	100 ROLLS
2	8 cm x 4 mtr	100 ROLLS
3	10 cm x 4 mtr	75 ROLLS
4	15 cm x 4 mtr	50 ROLLS

POP BANDAGE

DESCRIPTIONS:

Plaster bandages are fast and practical to use; soaked in water, the moistened bandage is applied to the limb, hardening rapidly, to form a resistant casting, which does not lose shape. Modeling plasters are used in Orthopedics as impression and modeling compounds



Usage Applications	Hospital & Clinical
Color	White
Feature	OrthoCare
Material	Imported POP Powdered
Sterility	Non Sterile
Brand	Osseo Care

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)


S.No	SIZE	MOQ
1	10 cm x 2.7 mtr	50 ROLLS
2	15 cm x 2.7 mtr	50 ROLLS

SOFT ROLLS (ORTHO ROLLS)



DESCRIPTIONS:

Generally used for primary wound dressings or is applied over a dressing, these bandage roll allow low adherent wound interface. These bandage roll have excellent absorbency, which promotes healing and also protects the wound.



Usage Applications	Hospital & Clinical
Color	White
Feature	OrthoDressing
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Quality	Premium

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	10 cm x 3 mtr	100 ROLLS
2	15 cm x 3 mtr	100 ROLLS

SURGICAL COTTON

DESCRIPTIONS:

Crepe bandage is suitable for use as support for Sprains and strains in joints and muscles. It also Works as a moderate compression bandage for Knee swelling, ankle swelling, and other relevant Injuries



Usage Applications	Hospital & Clinical
Color	White
Feature	Dressing
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Brand	Premium

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	25 GM	800 ROLLS
2	50 GM	400 ROLLS
3	100 GM	200 ROLLS
4	200 GM	100 ROLLS
5	400 GM	75 ROLLS

SURGICAL COTTON BALLS

DESCRIPTIONS:

For medical purposes, cotton pads are used to stop or prevent bleeding from minor punctures such as injections or venipuncture. They may be secured in place with tape. Cotton pads are also used in the application and the removal of makeup. Cotton pads are soft enough that they can be used to clean babies.



Usage Applications	Hospital & Clinical
Color	White
Feature	Dressing
Material	100% Cotton, Soft & Hygienic
Sterility	Non Sterile
Quality	Premium


AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

<i>S.No</i>	<i>SIZE</i>	<i>MOQ</i>
1	1 kg Cotton Balls	10 PKT

DISPOSABLE BED SHEETS WITH PILLOW COVER

DESCRIPTIONS:



Disposable Bed Sheets are made from soft non- woven material and crafted to provide protection against dust, allergies, and contamination. Ideal to be used in the clinic, hospitals, parlors, and spa, etc.

Usage Applications	Hospital & Clinical
Color	Navy Blue
Feature	Bed Care
Material	Non Woven
Sterility	Non Sterile
Quality	35 GSM

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	2.1 Mtr x 1.6 Mtr	100 PKT

UNDERPADS

DESCRIPTIONS:

Under pads, also known as Choux, are super absorbent pads designed to protect the mattress from urine damage. The pad is placed above or below the linens (personal preference) and absorbs leaking liquid.



Usage Applications	Hospital & Clinical
Color	Green
Feature	Pads
Material	100% Cotton
Sterility	Non Sterile
Quality	Premium

AVAILABLE SIZES :

(GST 12% & TRANSPORT EXTRA)

S.No	SIZE	MOQ
1	60 cm x 90 cm	50 PKT

NOTE:

**GST Extra Applicable (For Surgical Items
12%, For Rehabilitation Aids 5%)**

Transport Charges Extra

Payment Advance.

Thank you for your interest